


# FREE STANDING STERILIZERS

## vertical sliding doors

A wide range of free-standing sterilizers are available.

Chambers are built on the basis of three standard sections, **but we can quote other specific sizes upon request.**



Type	4307 S/D*	4707 S/D*	4710 S/D*	4713 S/D*	6710 S/D*	6713 S/D*	6719 S/D*
<b>DIMENSIONS OF THE CHAMBER</b>							
WIDTH (MM)	420	420	420	420	660	660	660
HEIGHT (MM)	350	700	700	700	700	700	700
LENGTH (MM)	620	620	1000	1300	1000	1300	1900
TOTAL VOLUME	92 liters	180 liters	287 liters	373 liters	455 liters	590 liters	870 liters
<b>STEAM GENERATOR</b>							
CAPACITY	24 liters	24 liters	24 liters	24 liters	24 liters	54 liters	67 liters
POWER	18 kW	18 kW	30 kW	30 kW	36 kW	36 kW	54 kW
<b>VACUUM PUMP</b>							
NOMINAL THROUGHPUT	55 m3/h	55 m3/h	55 m3/h	75 m3/h	75 m3/h	75 m3/h	106 m3/h
POWER	1,5 kW	1,5 kW	1,5 kW	2,2 kW	2,2 kW	2,2 kW	3,5 kW
<b>WEIGHT OF THE AUTOCLAVE</b>							
EMPTY	600 kg	700 kg	850 kg	950 kg	1150 kg	1400 kg	1650 kg
WITHH WATER	690 kg	880 kg	1137 kg	1323 kg	1605 kg	1990 kg	2520 kg
<b>TOTAL POWER</b>							
-	17 kW	22 kW	35 kW	35 kW	40 kW	40 kW	62 kW

\* S : simple door – D : double doors



Techniques bio médicales

Z.A. du Paradis, 2 av. Elie Baylac  
 95660 Champagne-sur-Oise - France  
 Tel. : +33 (0) 1 34 70 03 97 - Fax : +33 (0) 1 34 70 10 26  
 Internet : [www.sterichem.com](http://www.sterichem.com) - E-mail : [sales@sterichem.com](mailto:sales@sterichem.com)

